


# HeartTrends Quick User Guide

UG-C106 v1.4



## 1. Setup Phone (one-time)


- Download 'HeartTrends' App from Google Play Store or App Store
- Install the App
  - Accept each Permission with "Yes"
  - Android v8.0 or higher
  - Apple iOS v12.0 or higher
- Enter Username & Password (obtained from your Distributor)



"HeartTrends"

## 2. Wear (Wakeup) Sensor


- Adjust the strap on the skin around your bare chest (not the stomach).
- Wear the heart rate sensor by snapping on both ends of the belt
- Moisten the 2 rubberized belt electrodes for better connectivity & cleanliness (optional).
- Using an optical armband? See #8



Belt electrodes

## 3. Pair Sensor to App

- Press the **red** communication icon (upper right)
- Select your sensor from the list
- Press Connect
- Communication icon will turn **green**
- No need to pair this sensor using generic Bluetooth pairing



## 4. Take the Test

- Press "Start test"
- Enter Name and ID & optional medical information – important!
- Running Time Clock will appear
  - Option: Testing more subjects? Press **green** icon to disconnect & attach next device
- Wear the belt & sensor for **at least 20 minutes** (no need to remain connected to phone)
- Remove belt at end of test
- Return it to the technician



## 5. Retrieve Test Score

- Press the Retrieve Test button or **Red** Icon
- Wake up the sensor by **holding** the straps and connecting the snaps
- Select sensor from list. No response? Tap the sensor twice or unsnap & re-snap.
- If an old test resides on the sensor you will be prompted to submit it for analysis



## 6. Analyze your Test Score

- Physician should put score into context
- A score greater than 2.6 is high
- Scores less than 2.6 but greater than 2.0 is medium
- Scores below 2.0 is low. It may occur if subject was too active or stressed
- Suggest *repeating* the test in a quieter and more sedentary surrounding
- If 2<sup>nd</sup> test is also below 2.0, suggest consulting your physician



## 7. Review Test Results

- All test results may be viewed from the Review Test Results menu
- Detailed information obtained by pressing each test or the PDF icon
- Alternatively, log onto the HeartTrends website Clinic tab with your same login to view all test results in a tabular mode ([www.HeartTrends.com](http://www.HeartTrends.com))



## 8. Polar OH-1 armband: Manual Start/Stop required



SILVER BUTTON



LED STATUS INDICATOR

FUNCTION	SILVER BUTTON ACTION	LED FLASHING
ON (standby) / OFF	Long press	Once/second
RECORD start	Press twice (x2) quickly	Twice/second
RECORD stop	Press twice (x2) quickly	Once/second
Low battery	----	<b>Green + Red</b>